

GTK Optical Trocar

G T.K Disposable Optical Trocars are used for puncturing the abdominal wall, approaching the abdominal cavity, supplying gas and providing an access for surgical instruments in laparoscopic surgery.

CE 0197
ISO 13485 TÜVRheinland[®]
FDA 510K

Features:

Visible:

An endoscope can be inserted into an optical trocar to make the whole surgery process visible.

Reliability:

21 inspection steps on each trocar ensures high quality.

Safety:

The ergonomic design of the trocar thread not only makes instruments stable but also reduces harm to skin.

Design:

Pattern and color design is pleasing and comfortable to eyes.



Universality:

No need to have an adapter between 5mm and 15mm instruments. It is compatible to all kinds of surgical instruments.


Leak-proofness:

The universal seal can be replaced without removing the trocar cannula, ensuring 100% leak-proofness.

Detachability:

The sealing components can be detached to remove large tissue or organ.

What makes GTK Optical Trocars *stand out* ?



Product Name	Model	Product code	Number of Cannula	Diameter (mm)	Working Length(mm)	Qty Per Carton
Optical Trocars	CH-Optical	T05-75-CT-F	1 pc 5mm	5	75	20
		T12-75-CT-F	1 pc 12mm	12	75	20
		T05-100-CT-F	1 pc 5mm	5	100	20
		T10-100-CT-F	1 pc 10mm	10	100	20
		T12-100-CT-F	1 pc 12mm	12	100	20
		T12L-100-CT-F	1 pc 12mm	12	100	20
		T15-100-CT-F	1 pc 15mm	15	100	20
		T05-150-CT-F	1 pc 5mm	5	150	15
		T10-150-CT-F	1 pc 10mm	10	150	15
		T12-150-CT-F	1 pc 12mm	12	150	15
		2T052-75-CT-F	2 pcs 5mm	5	75	20
		2T052-CT-F	2 pcs 5mm	5	100	20
		2T102-CT-F	2 pcs 10mm	10	100	20
		2T122-CT-F	2 pcs 12mm	12	100	20
		2T122-75-CT-F	2 pcs 12mm	12	75	20
		2T12L2-CT-F	2 pcs 12mm	12	100	20
		3T52101-CT-F	2 pcs 5mm & 1 pc 10mm	5 & 10	100	10
		3T52121-CT-F	2 pcs 5mm & 1 pc 12mm	5 & 12	100	10
		3T5212L1-CT-F	2 pcs 5mm & 1 pc 12mm	5 & 12	100	10
		4T52102-CT-F	2 pcs 5mm & 2 pcs 10mm	5 & 10	100	10
		4T52122-CT-F	2 pcs 5mm & 2 pcs 12mm	5 & 12	100	10
		4T5212L2-CT-F	2 pcs 5mm & 2 pcs 12mm	5 & 12	100	10

GTK Bladed Trocars

G T.K Disposable Bladed Trocars are used for puncturing the abdominal wall, approaching the abdominal cavity, supplying gas and providing an access for surgical instruments in laparoscopic surgery.

CE 0197
ISO 13485 TÜVRheinland[®]
FDA 510K

Features:

Reliability:

21 inspection steps on each trocar ensures high quality.

Safety:

The ergonomic design of the trocar thread not only makes instruments stable but also reduces harm to skin.

Design:

Pattern and color design is pleasing and comfortable to eyes.



Universality:

No need to have an adapter between 5mm and 15mm instruments. It is compatible to all kinds of surgical instruments.

Leak-proofness:

The universal seal can be replaced without removing the trocar cannula, ensuring 100% leak-proofness.

Detachability:

The sealing components can be detached to remove large tissue or organ.

What makes GTK Bladed Trocars stand out ?



Product Name	Model	Product code	Number of Cannula	Diameter (mm)	Working Length(mm)	Qty Per Carton
Bladed Trocars	CH-A	T05-100-CT-A	1 pc 5mm	5	100	20
		T10-100-CT-A	1 pc 10mm	10	100	20
		T12-100-CT-A	1 pc 12mm	12	100	20
		2T052-CT-A	2 pcs 5mm	5	100	20
		2T102-CT-A	2 pcs 10mm	10	100	20
		2T122-CT-A	2 pcs 12mm	12	100	20
		4T52102-CT-A	2 pcs 5mm & 2 pcs 10mm	5 & 10	100	10
		4T52122-CT-A	2 pcs 5mm & 2 pcs 12mm	5 & 12	100	10

GTK Bladeless Trocars

G T.K Disposable Optical Trocars are used for puncturing the abdominal wall, approaching the abdominal cavity, supplying gas and providing an access for surgical instruments in laparoscopic surgery.

CE 0197
ISO 13485 TÜVRheinland®
FDA 510K

Features:

Visible:

An endoscope can be inserted into an optical trocar to make the whole surgery process visible.

Reliability:


21 inspection steps on each trocar ensures high quality.

Safety:

The ergonomic design of the trocar thread not only makes instruments stable but also reduces harm to skin.

Design:

Pattern and color design is pleasing and comfortable to eyes.



Universality:

No need to have an adapter between 5mm and 15mm instruments. It is compatible to all kinds of surgical instruments.

Leak-proofness:

The universal seal can be replaced without removing the trocar cannula, ensuring 100% leak-proofness.

Detachability:

The sealing components can be detached to remove large tissue or organ.

What makes GTK Bladeless Trocars stand out ?



Product Name	Model	Product code	Number of Cannula	Diameter (mm)	Working Length(mm)	Qty Per Carton
Bladeless Trocars	CH-F	T15-100-CT-F	1 pc 15mm	15	100	20 pcs
		T15-150-CT-F	1 pc 15mm	15	150	15pcs

What makes GTK Bladeless Pediatric Trocars stand out ?

GTK Disposable 3mm and 5mm trocars are used for puncturing the abdominal wall, approaching the abdominal cavity, supplying gas and providing an access for surgical instruments in laparoscopic pediatric surgery.

Features:

Reliability:

21 inspection steps on each trocar ensures high quality.

Safety:

The ergonomic design of trocar thread not only makes instruments stable but also reduces harm to skin.

Leak-proofness:

The universal seal can be replaced without removing trocar cannula, ensuring 100% leak-proofness.


Design:

Pattern and color design is amiable and comfortable to eyes.



What makes GTK Bladeless Pediatric Trocars stand out ?







The universal seal can be replaced without removing trocar cannula, ensuring 100% leak-proofness.

Rounded tip to ensure puncture safety.

The ergonomic design of trocar thread not only makes instruments stable but also reduces harm to skin.

Special obturator head design for better holding in puncture.

Product Name	Model	Product code	Diameter (mm)	Working Length(mm)	Qty Per Carton
Bladeless Trocars	CH-S	T03-55-CT-S	3	55	100
		T05-55-CT-S	5	55	100
		T05-100-CT-5	5	100	100