




G.T.K Medical



Veress Needle



Veress Needle

Safety

  Double warning system confirms proper insertion.

Quality

Sharp 304 stainless steel needle with strong polycarbonate handle.
16 quality inspection steps on every needle.

Convenience

Ergonomic handle for easy controlled insertion.



GT.K Medical Veress Needle/Insufflation needle is used in laparoscopic surgery for puncturing abdominal cavity and supplying CO₂ gas to create pneumoperitoneum.

Product Code	Working Length (mm)
VN-120	120
VN-150	150



CE 0197
ISO13485

©2013 Guangzhou T.K Medical Instrument Co., Ltd.

Add: A601 Guangzhou International Business Incubator, Guangzhou Science Park, Guangzhou, P.R. of China. Mobile: +86 139 0225 9487
Tel: +86 20 3229 0169 /3761 9648 Fax: +86 20 3229 0395/3761 6020 Http: //www.tk-medical.com E-mail: sales@tk-medical.com