



Wound Protectors/Retractors

Wound Protectors/Retractors are for disposable use in stabilizing and protecting wound incision in surgery.

Features and advantages:

360° wound incision expansion

Maximum exposure of a minimum incision;
Maximum surgical field;
Broaden surgical access;

360° wound incision protection

360° wound incision protection ;
Prevention of infection;
Minimize tissue damage;
Minimize bleeding of wound incision;

Time-saving

Save surgical time

Various models and diverse specifications

Three models--each model with diverse specifications for surgeon's choice.

Model	Specification	Diameter (mm)	Incision Size (mm)
Model A	WP-40	40	20~25
	WP-60	60	30~40
	WP-90	90	40~70
	WP-120	120	60~100
	WP-150	150	80~130
	WP-200	200	120~190
Model B	WP-40	40	20~25
	WP-60	60	30~40
	WP-90	90	40~70
	WP-120	120	60~100
	WP-150	150	80~130
	WP-180	180	90~170
	WP-200	200	120~190
	WP-230	230	140~220
	WP-250	250	150~240
	WP-280	280	170~260
	WP-330	330	220~310



Model B



Model B



G T.K Medical

WWW.TK-MEDICAL.COM

Model C



Wound Protectors/Retractors

Wound Protectors / Retractors



Model C



Model C



Model	Structure	Specification	ØD (mm)	ØD1 (mm)	ØD2 (mm)	L (mm)	Incision Size (mm)
Model C		WP-15C	15	30	36	20	10~20
		WP-25C	25	45	50	25	20~30
		WP-35C	35	70	70	25	30~40
		WP-60C	60	100	120	25	55~65
		WP-60TC	60	120	120	25	55~65
		WP-25SC	25	50	60	5	20~30
		WP-35SC	35	60	70	5	30~40
		WP-35LC	35	70	70	40	30~40
		WP-60LC	60	100	120	40	55~65
		WP-60LTC	60	120	120	40	55~65
		WP-60XLC	60	100	120	70	55~65

Quick Instructions for Use

Model A: Tear the Double-faced Adhesive Tape and adhere it onto the skin then it can be stabled and protected.

Model B: By rolling down the Rolling Ring the length of Protecting Film can be adjusted to fit for wound incisions of various lengths.

Model C: Fixed length and size ensures easier and convenient operation.

Statement: The products listed in this page are only sold in non-patented countries and regions.

©2023 Guangzhou T.K Medical Instrument Co., Ltd.

Add: No. 7, Jianye Six Road, East Zone of Guangzhou Economic and Technological Development District, Huangpu District, 510765, Guangzhou, China

Mobile: +86 139 2247 0081

+86 139 0225 9487

Tel: +86 20 3229 0169/3761 6021

Fax: +86 20 32290395

Http: //www.tk-medical.com

E-mail: sales@gtkmedical.com